

**Tues 4/26:**

- 1:39 - begin alignment procedure
- Laser height (above rail) = 19.51 +/- 0.1 cm
- SA is aligned to laser with retroreflection beside exit aperture
- 1:52 - Spectrum analyzer controller calibrated
  - Rise-time, offset, and amplitude following the instructions
  - Set to:
    - SWEEP EXPANSION: X1
    - AMPLITUDE: 250V
    - OFFSET: MID-RANGE
    - RISE-TIME: FULLY CCW (10 mS)
- 1:55 - Setting two channel oscilloscope:
  - CH 1:
    - Coupling: AC
    - BW Limit: OFF
    - Volts/Div: COARSE
    - Probe: 1X
    - Volts/Div: 1 V (Adjust via control in CH 1 panel)
  - CH 2:
    - Coupling: AC
    - BW Limit: OFF
    - Volts/Div: COARSE
    - Probe: 10X
    - Volts/Div: 50.0 V
- 2:12 - Adjust amplitude to see two peaks in oscilloscope and turn Volts/Div to 20 mV
- 2:38 - Change scope to 1ms/div and adjust amplitude and offset to have two clusters 8 divisions apart (using points that are in phase).
  - 0.125 FSR/div on scope
- 3:01 - Begin A) Polarization
  - Find degree of polarization
- 3:03 - Set up adjustable polaroid linear polarizer 70cm front of laser, with laser 80cm away
  - Wait for mode to change and then take measurements while changing angle of polarizer
- 4:08 - calculate degree of polarization w eq. 4:
  - $$\rho = \frac{I_{max} - I_{min}}{I_{max} + I_{min}} = 0.996 \pm 0.04$$
  - Need to graph  $P(\theta)$  vs.  $\theta$
- 4:12 - Begin B) Beam Diameter
- 4:14 - Exercise IV.8:
  - The direction of the beam is the z-axis
  - Write down the equation for the intensity  $I(x, y)$  of the beam as a function of the perpendicular distance

- $I(x, y) = \frac{P}{\pi r^2} = \frac{P}{\pi(x^2 + y^2)} = \frac{\epsilon_0 E^2}{2\pi(x^2 + y^2)}$
- 4:54 - Exercise IV.9:
  - Measure the beam profile with a power meter, and an adjustable knife edge
  - Mount the laser, power meter, and knife edge on the optical bench
  - For 1st measurement, place knife edge 1-2 cm from laser
  - Measure the total beam power  $P(x)$  (arbitrary units will do) as a function of the knife edge position,  $x$

#### Thurs 4/28:

- 1:30 - set up laser and SA
- Exercise IV.10:
  - $I(x, y) = c I_0 \cdot e^{-k(x^2 + y^2)}$
  - Where  $c$  and  $k$  are proportionality factors
  - we know that at  $x = y = 0$ , the intensity should be at its max, and that the intensity decays exponentially
- Power meter has uncertainty of  $\pm 1\%$  of total +  $\pm 5\%$  from detector

#### Tues 5/3:

- 1:40 - Set up laser and SA
- 1:41 - Exercise IV.1:
  - Using a python model, we got an equation to describe the data
    - beam diameter is the diameter at which the electric field has decayed to  $1/e$  of its maximum value, so that means that the intensity will have decayed to  $\frac{1}{e^2}$ .
    - Center at  $0.003229 \text{ cm}$
    - $\frac{1}{e^2} \cdot I_{max} = 114.9 \text{ uW}$
    - Position associated w this is at  $-0.015472 \text{ cm}$
    - $d_{beam} = 0.055 \text{ cm} + 0.015472 \text{ cm} = 0.070472 \text{ cm} \pm 0.0059 \text{ cm}$   
(8.3%)
- 2:40 - C) Beam Divergence
  - Wavelength =  $632.8 \text{ nm} = 6.328 \times 10^{-4} \text{ mm}$
  - $(178.4 \pm 1) - (30 \pm 0.4) = 148.4 \pm 0.25 \text{ cm} = l$
  - Trial 1:
    - Max (835) at  $1.01 \pm 0.001$
    - $1/e^2$  of max (113) at  $0.779 \pm 0.005$
  - $D_0 = d_{beam} = 0.070472 \pm 0.0059 \text{ cm}$
  - $D_1 = 0.231 \pm 0.019 \text{ cm}$  (8.2%)
  - Let  $2\Delta\theta = \theta = \text{divergence}$

- From wikipedia:  $\theta = 2\arctan\left(\frac{D_1 - D_0}{l}\right)$ 
  - From measurements:  $\theta = 1.08 \pm 0.18 \text{ mrad}$
- From manufacturer:  $\theta = 1.0 \text{ mrad}$
- From Fraunhofer diffraction:
  - $\theta = 2.44\left(\frac{\lambda}{D_0}\right) = 2.19 \pm 0.18 \text{ mrad}$
- 4:01 - D) Beam Waist:
  - Lens has  $f = 300 \pm 50 \text{ mm}$
  - Exercise IV.15:
    - $d = \left(\frac{4\lambda}{\pi}\right)\left(\frac{f}{D_0}\right) = 0.343 \text{ mm} \pm$
    - $L = 2\left(\frac{4\lambda}{\pi}\right)\left(\frac{f}{D_0}\right)^2 = 292.022 \text{ mm} \pm$

#### Thurs 5/5:

- 1:16 - set up laser
- D) Beam waist continued:
  - Set lens directly in front of laser
  - Place power meter 1 meter away from the lens
  - Take measurements of beam diameter in increments of 5cm, starting 5 cm in front of lens
- 4:42 - Having difficulty getting data that matches our expected result. The beam diameter constantly has a greater value that is not within uncertainty
  - Tried both sides of the lens to see if that would influence the data, but this did not seem to fix the problem
  - The intensity of the laser fluctuated from 760 uW - 815 u W, most likely due to the TEM mode transitions.
  - Beam diameter is affected by the maximum intensity of the laser, as we saw in part B, and finding the beam diameter consecutively over a course of time led to data that varied far too much
    - If the intensity increased while measuring, the beam diameter would increase and skew the data so that the beam waist could not be determined accurately.
    - The same would happen inversely if the intensity increased, but even with this the data collected was much larger than our theoretical value